



Ready to Use CHAPS Cell Extract Buffer

| Product Number | Amount/Size | Price |
|----------------|-------------|-------|
| #B0001 | 10 ml | \$25 |

Description

Ready to Use - Chaps Cell Extract Buffer is preferred for Western Blotting analysis using caspase antibodies.

Application

Total cell lysates are collected and lysed in CHAPS Cell extract buffer. Add protease inhibitors* to Chaps Cell Extract Buffer prior to use. After experimental treatment and washing, lyse cells in ice cold Chaps cell extract buffer. The lysate is processed by sonication and subsequent centrifugation at 12,500rpm to remove cell debris. The supernatant can be frozen at -20C until use for Western Blot analysis.

Storage Conditions

Store at -20°C. avoid repeated freeze/thaw cycle.

Companion Products

Caspase 3 Apoptosis Extracts (#E0003)
Caspase 8 Apoptosis Extracts (#E0008)
Caspase 9 Apoptosis Extracts (#E0009)
PARP Apoptosis Extracts (#S0001)
XIAP Apoptosis Extracts (#S0002)
BID Apoptosis Extracts (#S0003)

*Protease Inhibitors not included

For research purposes only. Not for use in humans.